



The Dolomite Centre Ltd.
Unit 1, Anglian Business Park, Orchard Road,
Royston, Hertfordshire, SG8 5TW, UK

T: +44 (0)1763 242491
F: +44 (0)1763 246125
E: dolomitesales@syrris.com
W: www.dolomite-centre.com

The Dolomite Centre Limited
Registered office: 27 Jarman Way,
Royston, Hertfordshire, SG8 5HW, UK
Company No. 05640084

INFORMATION SHEET

Part name	Miniature peristaltic pump (pack of 2)	Part number	3200024
-----------	--	-------------	---------

Description

The miniature peristaltic pump is a compact pump that produces flow at 0.45 ml/min. It has very low power consumption of 0.12W and is supplied with silicone tubing fitted.



Miniature peristaltic pump

Benefits

- Compact size
- Low power consumption
- Lightweight

	Specification	Value
1	Dimensions	30 x 12 x 14 mm
2	Flow Rate	0.45 ml/min
3	Discharge Pressure	0 – 50 kPa
4	Tube Size	OD 2.5 mm, ID 1.5 mm
5	Tube Material	Silicone
6	Motor	DC Geared Motor
7	Voltage	3V DC
8	Power Consumption	0.12 W
9	Weight	11 g

Geometry

